





- ▶ Automated needle retraction
- ▶ Effectively reduces the risk of needlestick injuries
- ▶ Quick and easy to use
- ▶ Activation can be achieved in a variety of ways
- ▶ May be activated before or after removal from patient
- ▶ Pre-removal activation reduces the risk of exposure to the contaminated sharp
- ▶ Activation of safety mechanism automatically covers syringe luer tip, reducing risk of exposure to medication, blood/body fluids
- ▶ Remains safe through disposal



(50 per box /400 per case)

NEEDLES ONLY				WITH 3mL SYRINGES			
Catalog #	Needle Size	Catalog #	Needle Size		Catalog #	Needle Size	
For syringes 3mL and smaller		For syringes 10mL and smaller					
82211	30G x 1/2"	85291	25G x 5/8"		82391	25G x 5/8"	
82091	25G x 5/8"	85211	25G x 1"		82311	25G x 1"	
82011	25G x 1"	85011	25G x 1 1/2"		85311	25G x 1 1/2"	
82031	23G x 1"	85231	23G x 1"		82331	23G x 1"	
		85021	22G x 1"		85321	22G x 1"	
		85031	22G x 1 1/2"		85331	22G x 1 1/2"	
		85041	21G x 1"		85341	21G x 1"	
		85051	21G x 1 1/2"		85351	21G x 1 1/2"	
		85061	20G x 1"		85361	20G x 1"	
		85071	20G x 1 1/2"		85371	20G x 1 1/2"	
		85081	18G x 1"		85381	18G x 1"	
		85091	18G x 1 1/2"		85391	18G x 1 1/2"	

#### OTHER AVAILABLE SAFETY PRODUCTS

- ▶ VanishPoint® Syringe
- ▶ VanishPoint® Blood Collection Tube Holder
- ▶ VanishPoint® Blood Collection Set
- ▶ VanishPoint® IV Catheter
- ▶ EasyPoint® Blood Collection Tube Holder with Needle



**RETRACTABLE**  
TECHNOLOGIES, INC.®

[rtisales@retractable.com](mailto:rtisales@retractable.com)

511 Lobo Lane • Little Elm  
Texas • 75068-5295 • USA  
Tel.: (972) 294-1770  
Fax: (972) 294-4400  
Toll Free: 1-888-703-1010

© 2015, 2025 Retractable Technologies, Inc. All rights reserved. M098 REV.9 EN, MAR 2025



**SAVES MONEY, TIME, AND LIVES**



[www.retractable.com](http://www.retractable.com)




**RETRACTABLE**  
TECHNOLOGIES, INC.®

# easypoint™

RETRACTABLE NEEDLE

**SAVES MONEY, TIME, AND LIVES**




## easy to use

- Minimal change in injection technique
- Allows for single-handed activation
- Clinician's hands and fingers remain behind the exposed sharp at all times



1 Securely attach EasyPoint needle to a syringe. Administer medication according to institutional policy.



2 Apply pressure to the color-coded tab and the safety chamber using a squeezing motion, allowing the needle to automatically retract into the safety chamber.



3 Dispose in an appropriate sharps container.